

User Guide for Training Classification Models

```
Starting training for 500 epochs...

Epoch  GPU_mem  box_loss  seg_loss  cls_loss  dfl_loss  Instances  Size  57/57 [00:11<00:00, 5.04it/s]
1/500   3.75G    3.068    6.457    4.865    4.122     23         640: 100%|  9/9 [00:02<00:00, 3.08it/s]
      Class  Images  Instances  Box(P  R  mAP50  mAP50-95)  Mask(P  R  mAP50  mAP50-95)
      all    337     349    0.000462  0.285  0.0128  0.00522  0.00017  0.126  0.00542  0.00131

Epoch  GPU_mem  box_loss  seg_loss  cls_loss  dfl_loss  Instances  Size  57/57 [00:08<00:00, 6.33it/s]
2/500   3.74G    2.894    3.341    4.223    3.526     16         640: 100%|  9/9 [00:02<00:00, 3.89it/s]
      Class  Images  Instances  Box(P  R  mAP50  mAP50-95)  Mask(P  R  mAP50  mAP50-95)
      all    337     349    0.381    0.0197  0.0146  0.00583  0.359  0.0263  0.0137  0.00884

Epoch  GPU_mem  box_loss  seg_loss  cls_loss  dfl_loss  Instances  Size  57/57 [00:08<00:00, 6.57it/s]
3/500   3.79G    2.572    2.963    3.696    3.038     26         640: 100%|  9/9 [00:02<00:00, 3.70it/s]
      Class  Images  Instances  Box(P  R  mAP50  mAP50-95)  Mask(P  R  mAP50  mAP50-95)
      all    337     349    0.105    0.173  0.046  0.0177  0.115  0.147  0.0452  0.0223

Epoch  GPU_mem  box_loss  seg_loss  cls_loss  dfl_loss  Instances  Size  57/57 [00:09<00:00, 6.32it/s]
4/500   3.79G    2.242    2.645    3.118    2.766     12         640: 100%|  9/9 [00:02<00:00, 4.06it/s]
      Class  Images  Instances  Box(P  R  mAP50  mAP50-95)  Mask(P  R  mAP50  mAP50-95)
      all    337     349    0.195    0.293  0.155  0.0717  0.189  0.289  0.156  0.0817
```

Note:

The following example shows how the model training was done and how to recreate it. To save time and storage space, the batch size and epochs were reduced, and only a small fraction of the training data was included.

Step 1:

Clone the project repo:

git clone https://github.com/tchiang0/data_515_brain_tumor_computer_vision.git

Step 2:

Install pip if not already installed:

```
python -m pip install
```

Step 3:

Run requirements.txt to ensure all dependencies exist:

```
pip install -r requirements.txt
```

Step 4:

Navigate into the data_515_brain_tumor_computer_vision project directory and then into data_515_brain_tumor_computer_vision folder:

```
cd data_515_brain_tumor_computer_vision/data_515_brain_tumor_computer_vision
```

Step 5:

Run the following command from inside the data_515_brain_tumor_computer_vision folder to train the two classification models (one to determine if the uploaded image is a brain scan, one to determine if the uploaded brain scan has a tumor).

```
python model_training/model_building.py
```


```
PS C:\Users\Elliott\Documents\GitHub\data_515_brain_tumor_computer_vision> cd .\data_515_brain_tumor_computer_vision\
PS C:\Users\Elliott\Documents\GitHub\data_515_brain_tumor_computer_vision\data_515_brain_tumor_computer_vision> python model_training/model_building.py
Downloading https://github.com/ultralytics/assets/releases/download/v8.1.0/yolov8n-cls.pt to 'model_training\yolov8n-cls.pt'...
100% | 5.30M/5.30M [00:00<00:00, 22.6MB/s]
New https://pypi.org/project/ultralytics/8.1.27 available 🟡 Update with 'pip install -U ultralytics'
Ultralytics YOLOv8.1.10 🚀 Python-3.9.10 torch-2.2.0+cu121 CUDA:0 (NVIDIA GeForce RTX 3080, 10240MiB)
engine\trainer: task=classify, mode=train, model=model_training\yolov8n-cls.pt, data=model_training/sample_data/is_scan_data, epochs=1, time=None, patience=50, batch=3, imgsz=640, save=True, save_period=-1, cache=False, device=None, workers=8, project=None, name=is_brain_scan, exist_ok=False, pretrained=True, optimizer=auto, verbose=True, seed=0, deterministic=True, single_cls=False, rect=False, cos_lr=False, close_mos
```

Step 6:

Navigate to the folder where the weights and results are stored:

runs/classify/is_tumor

runs/classify/is_brain_scan



Final notes:

These steps are also executed automatically when running program tests. To see which training data is used, navigate to the model_training/sample_data folder.